

# NCA Tutorial

Yingbo Ma

30nd November 2018

## 1 Introduction

This is an introduction to `NCA.jl`, a software for noncompartmental analysis (NCA). In this tutorial we will show how to use `NCA.jl` to analysis data.

### 1.1 Installation

Currently, `NCA.jl` is a submodule in `PuMaS.jl`, so you only need to install `PuMaS.jl`, and everything will be ready to go.

### 1.2 Getting Started

To load the package, use

```
using PuMaS.NCA
```

First, let's load the example NCA data inside `PuMaS.jl`. This data have 24 individuals, and each of them has 16 data points. Here, we are creating two functions `concs(i)` and `times(i)` which give the concentration timeseries of individual `i`.

```
using PuMaS, CSV
```

```
const example_nca_data = CSV.read(example_nmtran_data("nca_test_data/dapa_IV"))
concs(i) = Float64.(example_nca_data[:, :CObs])[16(i-1)+1:16*i]
times(i) = Float64.(example_nca_data[:, :TIME])[16(i-1)+1:16*i]
```

## 2 Efficient Computation of Multiple NCA Diagnostics

### 2.1 AUC and AUMC

We can compute the area under the curve (AUC) from the first observation time to infinity

```
auc(concs(1), times(1))
```

```
263.792662196049
```

```
auc(concs(1), times(1), method=:linuplogdown)
```

```
257.8586273987722
```

the keyword argument `method` can be `:linear`, `:linuplogdown`, or `:linlog`, and it defaults to `:linear`. This is a simple interface, however it is not efficient if you want to compute many quantities. The recommended way is to create an `NCAdata` object first and then call the respective NCA diagnostic on the data object:

```
nca = NCAdata(concs(1), times(1))
auc(nca)
```

```
263.792662196049
```

Let's create `NCAdata(concs(i), times(i))` for each individual `i` from 1 to 24. To do that, we would broadcast the whole expression over the individual indices we'd want to use, like:

```
ncas = @. NCAdata(concs(1:24), times(1:24));
```

`AUClast` is the area under the curve from the first observation to the last observation. To compute `AUClast` on the first individual, one would do:

```
ncas = @. NCAdata(concs(1:24), times(1:24), dose=5000.)
auc(ncas[1], auctype=:AUClast)
```

```
(246.932184, 0.0493864368)
```

Note that if `dose` is provided, `auc` will return a tuple, which is in the form of (AUC, normalized AUC).

Or to compute the AUC on every individual, one would do:

```
@. auc(ncas, auctype=:AUClast)
```

```
24-element Array{Tuple{Float64,Float64},1}:
 (246.932184, 0.0493864368)
 (302.24594, 0.060449188)
 (288.579555, 0.05771591100000001)
 (333.80397999999997, 0.066760796)
 (129.06118400000003, 0.025812236800000006)
```

```
(291.9510375, 0.0583902075)
(333.99432249999995, 0.06679886449999999)
(259.9668375, 0.0519933675)
(233.64285700000002, 0.0467285714)
(242.71946250000002, 0.048543892500000005)
⋮
(196.6721735, 0.039334434700000004)
(319.29748499999994, 0.06385949699999999)
(185.39912850000002, 0.037079825700000006)
(403.216205, 0.080643241)
(202.639537, 0.0405279074)
(222.76995, 0.04455399)
(364.75571750000006, 0.07295114350000001)
(265.6632825, 0.05313265649999999)
(156.80645950000002, 0.031361291900000005)
```

One can also compute AUC on a certain interval. To compute AUC on the interval  $[10, \infty]$  on the first individual

```
auc(ncas[1], interval=(10,Inf));
```

In many cases, the AUC commands may need to extrapolate in order to cover the desired interval. To see the percentage of extrapolation ( $\frac{\text{extrapolated AUC}}{\text{Total AUC}} \cdot 100$ ), you can use the command:

```
auc_extrap_percent(ncas[1])
```

```
6.391564517256502
```

Area under the first moment of the concentration (AUMC) is

$$\int_{t_0}^{t_1} t \cdot \text{concentration}(t) dt. \quad (1)$$

The interface of computing AUMC is exactly the same with AUC, and one needs to change `auc` to `aumc` for calculating AUMC or related quantities. For instance,

```
aumc(ncas[1])
aumc_extrap_percent(ncas[1])
```

```
59.47554145363463
```

## 2.2 Terminal Rate Constant ( $\Lambda Z$ )

The negative slope for concentration vs time in log-linear scale is the terminal rate constant, often denoted by  $\Lambda Z$ . To compute  $\Lambda Z$ , one can call

```
lambdaz(ncas[1])
```

```
(lambdaz = 0.03876710923615265, points = 5, r2 = 0.7444781209375907)
```

`lambdaz` returns a tuple in the form of ( $\Lambda Z$ , the number of data points used,  $r^2$ ), where  $r^2$  is the coefficient of determination from the linear regression of the slope determination. By default, it checks last 10 or less data points, one can change it by providing the keyword `threshold`, e.g.

```
lambdaz(ncas[1], threshold=15)
```

```
(lambdaz = 0.03876710923615265, points = 5, r2 = 0.7444781209375907)
```

One can also specify the exact data points by passing their indices

```
lambdaz(ncas[1], idx=[10, 15, 16])
```

```
(lambdaz = 0.03876710923615265, points = 5, r2 = 0.7444781209375907)
```

## 2.3 Simple functions

`T_max` is the time point at which the maximum concentration (`C_max`) is observed, and they can be computed by:

```
tmax(ncas[1])  
cmax(ncas[1])  
cmax(ncas[1], interval=(20, 24))
```

```
(0.653632, 0.0001307264)
```

note that `cmax` returns `C_max` and normalized `C_max` if `dose` is provided.

`T_last` is the time of the last observed concentration value above the lower limit of quantization (LLQ), and the corresponding concentration value is (`C_last`). They can be computed by the command

```
tlast(ncas[1])  
clast(ncas[1])
```

```
0.653632
```

By default, LLQ is 0, and concentrations that are below LLQ (BLQ) are dropped. If you wish to keep BLQ data and change LLQ to 0.5, you can do

```
NCAdata(concs(1), times(1), llq=0.5, condblq=:keep)
```

```
NCAdata:
  concentration: Float64
  time:         Float64
  auc:          Float64
  aumc:         Float64
  λz:           Float64
  dose:         Nothing
```

The half-time can be computed by:

```
thalf(ncas[1])
```

```
17.879774742490834
```

One may need to interpolate or to extrapolate the concentration-time data. For example, if you wanted to interpolate the concentration at  $t = 12$  using linear interpolation, you would do:

```
NCA.interpextrapconc(ncas[1], 12., interpmethod=:linear)
```

```
0.911367
```

interpmethod can be :linear, :linuplogdown, or :linlog.